



Subject: Performance of Beijing Lepu Medical Technology Co., Ltd's SARS-CoV-2 Antigen Rapid Test Kit are theoretically not be affected by recently discovered variants including the United Kingdom variant and the South African variant.

Dear valued customers:

We, Beijing Lepu Medical Technology Co., Ltd, would like to inform you that our SARS-CoV-2 Antigen Rapid Test Kit remain suitable for the detection of SARS-CoV-2 antigen even in the emergence of newly discovered variants including the United Kingdom variant and the South African variant.

According to our investigation, several site mutations have occurred of the spike protein at positions of N501Y, E484K, K417N/T for SA mutant strain B.1.351(501Y.V2), and of N501Y, P681H, 69-70del、E484K for UK mutant strain b.1.1.7(501Y.V1), (From Guangdong Provincial Center for Disease Control and Prevention) and of N501Y、E484K、K417N/T for Brazil mutant strain P.1(501Y.V3). Since the recognition site of the raw materials used in our antigen test are the nucleocapsid protein different from mutation sites, we expect our products are theoretically able to detect variants including UK, Brazil, and SA.

Meanwhile, we will promptly communicate any updates regarding SARS-CoV-2 Antigen Rapid Test Kit. In addition, we will continue our efforts to comply with high quality management standards and to maintains a consistent high quality management system to ensure customer's satisfaction and product safety. If you have any questions, please contact our sales representative.

Beijing Lepu Medical Technology Co., Ltd  
2021-02-19



北京乐普医疗科技有限责任公司

BEIJING LEPU MEDICAL TECHNOLOGY CO.,LTD.

注册地址：北京市昌平区科技园区超前路37号  
生产地址：北京市昌平区科技园区超前路37号7-1楼  
售后服务：北京市昌平区科技园区超前路37号7-1楼  
电 话：010-80123111

传 真：010-80123100  
邮 编：102200  
网 址：www.leputechnology.com  
咨询电话：400 060 1160